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Operation

Application

The centrifuge is only designed for centrifuging substances and substance mixtures with a density of max. 1.2 g/ml. Through the generation of centrifugal force, they separate substance mixtures or change their mixture ratio.

Safety instructions

- The centrifuge is designed to the state of the art and is operationally reliable; however, it can constitute a source of danger for the user or a third party if it is used improperly by untrained personnel or for purposes for which it was not designed.
- Before starting the centrifuge, read the operating instructions carefully.
- In addition to an operating and service manual and the binding regulations for accident prevention, the recognized industrial regulations for safety and professional work must be observed.
- The operating and service manual must be supplemented by accident prevention and environmental protection instructions based on national regulations prevailing in the country of the user.
- Place the centrifuge on a sturdy base.
- Load the centrifuge evenly.
- The maximum speeds printed on rotors must not be exceeded.
- Only use accessories for the centrifuge approved by the manufacturer.
- Centrifuge vessels may only be filled to the maximum level stated by their manufacturer.
- Do not operate the centrifuge in explosion-hazard rooms.
- If the centrifuge is used for hazardous substances and substance mixtures as listed in the accident prevention regulations, the operator shall take suitable steps to remain within the regulations.
- It is not permitted to operate this centrifuge using very corrosive substances which could have a detrimental effect on mechanical rotor strength.

- Should corrosion occur as a result of lack of care, the manufacturer will not accept any warranty claims.
- Any centrifugal rotors showing evidence of significant corrosion or even mechanical damage may not be used.
- Repairs may only be performed by persons authorized by the manufacturer.
- Before opening the housing, unplug the power cord from the mains socket.
- Only use original spare parts and permitted original accessories.
- The centrifuge is a Group 3 device in compliance with the Medizinische Geräteverordnung (German Medical Device Ordinance). It complies with the "Gesetz über technische Arbeitsmittel" (Law concerning technical working devices) dated 24 June 1968 with the supplement of 13 August 1979.

Commissioning

- Place centrifuge on a firm surface in a suitable location.
- Check whether the mains voltage is identical with the information on the nameplate.
- Plug the power cord to the mains socket.
- Turn mains switch into position "I", rotation indicator (function check) illuminates.
- After rotation indicator has gone off, open cover by pulling up safety clamp.

This cover can be opened only, if the centrifuge is switched on and the rotor at a standstill.
If this is not possible, see section on "Emergency Release".

Before first operation of this centrifuge, remove transportation supports (see instruction sheet). Keep transportation supports and original packaging in a safe place.

Fitting a Rotor (see Changing a Rotor)

Load Centrifuge

Only rotor and accessory combinations approved by the manufacturer for this centrifuge may be used (see "Rotors and Accessories" at the end of these operating instructions)!

During centrifugation all rotor compartments must be filled. Mixed rotor charging is permissible, if opposite compartments are filled with identical suspended trays.

- Free-swinging suspended trays must be attached to rotor journals such that they are easily swung (carry out swing test).
- Insert centrifugation vessels - which have been evenly filled using best visual judgement - at opposite locations into the rotor to avoid any imbalance.
- Close cover and push safety clamp forward. This safety clamp must audibly snap home at the two collar screws in the top of the casing. Control lamp will illuminate for three seconds.

 Centrifuge can only be started when cover is closed.

Operating Centrifuge for Preselected Time Periods

- Set timer to centrifuging time required.

 For short centrifuging times, first set timer to a longer centrifuging time to wind up the clock. Then set to actual centrifuging time required.

- Press start key, rotation indicator illuminates.
- Using speed selector and speed display, set maximum speed required (The figures shown on this speed selector are for guidance only).

 When carrying out centrifuging runs in series and at the same speed, only centrifuging times must be reset.

When the centrifuge has been loaded unevenly, its drive will cut out immediately. The rotation indicator will flash every 0.5 seconds. → When rotor has come to a stop, open cover and load rotor evenly.

Impulse Operation

- Set timer to "0".
- Press start key and hold for as long as this centrifuging run is to continue. Rotor will turn at the speed preselected via its speed selector.

Stopping the Centrifuge

- At the end of a preset centrifuging time, the timer will shut down the centrifuge. The rotor will be braked and slowed down until it stops.

 The centrifugation process can be interrupted at any time by resetting timer or speed selector to "0".

- When rotor has come to a stop, open cover and take out centrifugation material.
- Set mains switch to position "0".

Changing a Rotor

Removal

- Open cover
- Undo rotor tightening nut, using spanner wrench, and turn until lift-off pressure point is reached. When this lift-off pressure point has been overcome, the rotor will come off its motorshaft cone. Turn tightening nut until rotor can be lifted off its motor shaft.
- Lift rotor off its motor shaft.

Fitting

- Select rotor and mount on motor shaft. (Rotors are labelled with their respective nominal speed).

The motor shaft driver pin must snap into the groove at the bottom of the rotor.

- Tighten nut.

Check rotor for secure mounting.

Emergency Release

- If a power cut or any other type of fault occurs on the centrifuge whilst in operation, its locking system will remain blocked.

Open cover only when centrifuge has stopped!

- Pull up safety clamp.
- For emergency release purposes, there is a bore in the top front end.
 - Insert release pin E003 (supplied) into this bore until it reaches its stop position.
 - Press cover downwards and simultaneously move release pin to the right, until the locking mechanism is released.
 - Open cover

Centrifugation of Substances with a Density higher than 1.2 g/ml

Rotors are designed such that, at their respective rated maximum speed, they are capable of centrifuging substances with an average homogenous density of 1.2 g/ml maximum.

Substances with a higher density can only be centrifuged at reduced speeds.

Permitted speeds can be calculated using the following formula:

$$\text{Reduced speed (nred)} = \sqrt{\frac{1,2}{\text{higher density}}} \times \text{maximum speed}$$

$$\text{e. g.: nred} = \sqrt{\frac{1,2}{1,6}} \times 4000 = 3464 \text{ min}^{-1}$$

If you have any queries, ask the manufacturer for advice.

Calculating Relative Centrifugal Acceleration (RCA)

The relative centrifugal acceleration can be calculated as follows:

$$\text{RCA} = \left(\frac{n}{1000} \right)^2 \times r \times 1.118$$

$n = \text{min}^{-1}$ (speed)
 $r = \text{mm}$ (radius)

For a radius, use the distance between the journal shaft center line and the bottom of the centrifugal vessel (see "Rotors and Accessories" at the end of these operating instructions).

Care

- Centrifuge, rotor and sleeves must be cleaned regularly, using a weak solution of soapy water if necessary.
- In the event of **glass breakage** remove fragments and leaked centrifuging material from centrifugal chamber or from rotor.
For this task, the centrifugal chamber can be taken out of the centrifuge.

Rubber inserts in the suspended trays must be replaced as glass fragments embedded in these rubber inserts cause further glass breakage!

- If infectious centrifuging material enters into the centrifugal chamber or rotor, these are to be disinfected after cleaning.

Care of Rotors and Suspended trays (weekly)

To protect rotors and suspended trays against corrosion, they are to be cleaned using a weak solution of soapy water, if required.

For additional corrosion protection, aluminium rotors or suspended trays are to be smeared lightly with a suitable anti-corrosion agent.

Suspended trays showing evidence of cracks or deformations must be immediately replaced in complete sets!

Grease Journals (if required)

To keep journals on free-swinging rotors in gliding condition, thus achieving an even swing-out of inserted suspended trays, these journals must be greased with a suitable lubricating agent (Hettich 4051), if required.

Non-greased journals may cause operating faults (imbalance).

Check Sealing Ring

The sealing ring along the top edge of the centrifugal chamber is to be checked for damage (e.g. cuts, cracks, etc) and replaced if necessary

User Check (monthly)

Users need to check and ensure that the centrifuge is not visibly damaged. In particular, the following points need to be checked:

- motor mounting
- true running of motor shaft
- fault-free condition of rotors and suspended trays (no deformations or malformations)
- bolting

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Please enter the following information in the spaces provided:

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- Inventory no.
- Inspection no.
- Location of appliance

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Scope of Delivery

The following accessories will be included with this centrifuge:

	Sales no.
1 Connection cable (230 V)	4718
(115 V)	6021
2 Carbon brushes	E854
2 Fuse inserts (230 V)	
(F2,5 A) (5 x 20)	2386
or	
2 Fuse inserts (115 V)	
(T3,5 A) (5 x 20) UL/CSA	E770
1 Fuse insert	
(T 500 mA) (5 x 20)	E864
1 Spanner wrench	6327
1 Release pin	E003
1 Operating instruction	0034
1 Transportation support	
instructions	110.26.06.02

Rotor(s) and their relevant accessories will be custom supplied depending on requirements as stated in the order.

Maintenance/Repair

Repairs may only be carried out by persons authorised by the manufacturer.

Disconnect the centrifuge from the mains supply before opening the housing!

Following maintenance work, a function check is to be carried out.

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Change Fuses

Unplug the power cord from the mains socket.

The fuses are located inside the connector.

- Fold down connector flap.
- Press spring on fuse holder towards this fuse holder and withdraw holder complete with fuse.
- Change fuse.

Checking/Changing Carbon Brushes

The life of the carbon brushes mainly depends on the frequency of use of the centrifuge.

Check the carbon brushes at least every six months. In order to avoid damage to the motor or performance losses, always change the carbon brushes in pairs when they have worn down to 3 mm.

- Open cover.
- Unplug the power cord from the mains socket.
- Remove rotor (see Changing a Rotor)
- Take out centrifugation chamber from centrifuge.
- Disconnect earth conductor from centrifugation chamber.
- Remove the screw caps on the side of motor and remove carbon brushes.
- Replace carbon brushes (min. length of 3 mm) with new ones.

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- Reassemble centrifuge following above instructions in reverse order.

When replacing the rotor, make sure that the driving pin of the motor shaft engages in the rotor slot.

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- After repair, check whether:
 - all the plug-on connections have been replaced
 - all the screws are tight
 - all parts were reassembled
- Perform functional test.

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Changing a Motor or Changing Motor Mounting

- Open cover.
- Unplug the power cord from the mains socket.
- Remove rotor (see Changing a rotor)
- Remove centrifugation chamber in an upward direction.
- Disconnect earth conductors from centrifugation chamber and motor.
- Disconnect motor supply cable from electronic assembly.
- Undo hexagonal nut located at bottom of casing.
- Check motor mounting (rubber/metal mounting) for wear and, if necessary, replace as a complete set.
- Replace motor.
- Reassemble centrifuge following above instructions in reverse order.

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Mount motor such that it is centrally located inside centrifugation chamber.

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Changing Control Panel Elements

- Open cover.
- Unplug the power cord from the mains socket.
- Remove rotor (see Changing a rotor).
- Remove centrifugation chamber in an upward direction.
- Disconnect earth conductor from centrifugation chamber.
- Electrically disconnect relevant elements.
- Undo element at front or rear of control panel.
- Remove element from front or rear of control panel.
- Replace element.
- Reassemble centrifuge following above instructions in reverse order.
- Perform functional test.

Fault in Speedometer System

If during rotation a fault occurs in the speedometer system or an electronic assembly, the indicator light will start to flash at 3 second intervals.

The drive unit will shut down, the motor will run out without braking.

After 3 minutes, it is possible to open the cover.

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Open cover; set mains switch to "O"/"I". Indicator light will go off.

If the fault cannot be remedied in this way, contact your local customer service agent.

Functional Test

After repair, check whether:

- all the plug-on connections have been replaced
- all the screws are tight
- all parts were reassembled

After every repair, check the centrifuge for compliance with the characteristics listed in the table below.

The test readings may only deviate slightly from the characteristics in the table.

	Characteristics	
Rated voltage	230 V 50/60 Hz	115 V/60 Hz
Earth conductor resistor	< 0,3 Ω	< 0,3 Ω
Insulation resistor	> 2 M Ω	> 2 M Ω
Leak current	< 0,75 mA	< 0,5 mA

Finally, a function check is to be carried out (cf. "Technical Data").

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Customer Service


If any faults or disturbances which you cannot remedy yourself occur, please contact your nearest Hettich customer service centre or your nearest agents.

Transport

If you send a defective unit for repair to the manufacturers, the customer service centre or an agent, first re-install transportation supports. When shipping this centrifuge, use its original packaging only.

Schaltschema
UNIVERSAL 16A 230 V

Wiring Diagram
UNIVERSAL 16A 230 V

P1	Feldstrom (Bremsel) field current (brakes)
P2	Drehzahlanzeige speed indicator
P3	Drehzahl-nas. max. speed
P4	Bremschopper Ref. ref. value of brake chopper

X3, X4	S13, S14	M1	Rotar rotor			S15	Unwuchtschalter imbalance switch
	S12, S13		Netz-Sicherung F2,5A mains fuses F2,5A			S1	Netz EIN/AUS mains ON/OFF
	X9		Gerätekicker connector plug			S17	Sicherung T500mA fuse T500mA
	Z1		Funkentstörfilter radio suppression filter	X8	S18	S3	Start-/Kurzlaufschleife start/short cycle key
X2	S12	Y1	Deckelschalter/ Entriegelungsmagnet lid locking/ unlatching magnet			H1	Rotationsanzeige indicator/rotor running
						S4	Zeitschaltuhr line switch
X6	S16	R1	Drehzahlpoti potentiometer/speed	X7	S17	P5	Drehzahlanzeige speed indicator
Leistungs- teil A1 power amp. A1	Stecker- plug	Netzsch- schalter switch	Bezeichnung designation	Leistungs- teil A1 power amp. A1	Stecker- plug	Netzsch- schalter switch	Bezeichnung designation

rotor for setting imbalance values: rotor 1323 with buckets No.1301
 imbalance: run $\leq 10g$
 break $\leq 15g$

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